FT18 REVERSE OSMOSIS/ULTRAFILTRATION UNIT



The Armfield RO/UF unit is the smallest known industrial membrane system of its type in the world. It is designed specifically for laboratory scale testing of both reverse osmosis and ultrafiltration techniques.

features

- Dual operation, RO/UF
- Uses proprietary membranes available from PCI
- Clean-in-place facility
- Mobile and self-contained operation
- Computer linked for data logging and analysis

benefits

- Rapid changeover from R0 to UF and back
- Easily cleaned and maintained









Miniature-scale R & D technology

Head Office:

Armfield Limited Bridge House, West Street, Ringwood, Hampshire. BH24 1DY England

Telephone: +44 1425 478781 Fax: +44 1425 470916 E-mail: sales@armfield.co.uk

U.S. Office:

Armfield Inc.
436 West Commodore Blvd (#2)
Jackson, NJ 08527
Telephone: (732) 928 3332
Fax: (732) 928 3542
E-mail: info@armfieldinc.com



For more product information please visit: www.explorearmfield.com







twitter.com/Armfield_IFT facebook.com/Armfielduk linkedin.com/companies/armfield-limited youtube.com/user/armfieldUK explorearmfield.wordpress.com

Description

The FT18 is fully self-contained in a mobile cabinet. Comprehensive instrumentation is provided with digital meters displaying temperatures and system pressure. In addition to the membrane module, the cabinet also houses feed and permeate tanks, a plate heat exchanger and positive displacement pumps. Connection to cooling water and a suitable electrical supply is all that is required for normal operation.

All parts in contact with the process fluid are constructed from hygienic design materials such as stainless steel. PTFE and silicone rubber.

The RO/UF unit has an efficient clean-in-place facility and is easily dismantled for servicing or maintenance.

It uses a PCI MICRO 240 membrane module which accommodates all PCI membrane types. The unit is supplied with a pack of both R.O. membranes (AFC99) and U.F. membranes (ES625)

(other membranes not supplied by Armfield).

Membranes can be installed or removed rapidly for use in either process.

The filtration membranes differ only in length and quantity from those installed in the full-scale plant. Using small quantities of process fluid, operating parameters can be established and then optimised to effect a confident transition to full-scale.

Concentration/clarification/fractionation and standardisation of: milk whey fruit juices vegetable juices alcoholic beverages

Enzymes, antibiotics and organic acids



Removing the membrane from the module

Treatment of effluent

*PCI invites contact regarding membrane information and can provide details of local suppliers: Tel: +44 256 896966 Fax: +44 256 893835

Specifications

Flowrates: 30, 18, 15 and 9lpm Pressures: Variable up to 55 bar PCI MICRO-240 Membrane Module suitable for two

PCI* membranes

Diameter: 1.25cm
Area: 0.024m²
Tubeside volume: 75ml
Typical permeate flow: 5-50ml/min
Feed tank capacity: 15 litres
Permeate tank capacity: 15 litres

Reauirements

Electrical Supply: Three Phase (see Ordering Codes)

Water supply: 10 lpm @ 12°C: a normal labora-

tory cold water supply will suffice.

Ordering Codes

FT18-C: 415V/3ph/50Hz FT18-D: 208V/3ph/60Hz FT18-E: 380V/3ph/50Hz FT18-F: 220V/3ph/60Hz

(all units 2.2kW)

Optional Accessories

FT18-90IFD-USB: Data Logging Accessory

Requires user supplied PC with USB interface running Windows XP or later.

FT18-15: Splash screen

Shipping specification

Volume: 1.4 m³ Gross weight: 220kg

Overall dimensions

Height: 1.1m Width: 0.7m Depth: 0.7m

The Armfield range includes HTST/UHT/aseptic systems, carbonator/filler/cappers, spray dryers/chillers, multifunction batch processors, ice cream freezers, margarine crystallisers, extractors, edible oils processors and more. For further information about our products and services, or to book a trial at one of our trials facilities, please contact us.