



BUTTER CHURN AND BUTTER MAKING ACCESSORIES



FT21/21A Provisional

The Armfield Butter Churn is designed to provide practical training in the preparation of butter by the traditional method.

PROCESS CAPABILITIES

- 'phase inversion' of the oil/water emulsion which occurs in churning
- investigation of process temperature and agitation time in the agitation process
- production of butter by various methods for quality and analytical determinations
- mass balance quantative yield of butter from various types of raw milk



DESCRIPTION

The metal churn paddle is electrically driven. A multi-speed electric motor is mounted on the side of the vessel, which drives the paddle agitator. The paddle may be removed for cleaning. A batch of up to 5 litres of cream can allow production of more than 1kg of butter, and all of the important stages of this 'colloid inversion' process can be adequately observed.

TECHNICAL SPECIFICATIONS

Total churn volume: 5 litres Batch size: 5 litres 4 litres

5 litres (max) 4 litres (suggested) 1 litre (min)

SERVICES REQUIRED

Electrical supply: FT21-A: 220-240V/1ph/50Hz *FT21-B:* 120V/1ph/60hZ

SHIPPING SPECIFICATIONS

Volume: Gross weight:

0.03m³ 10kg

BUTTER MAKING ACCESSORIES - FT21A

The accessories have been selected to enable a finished product to be produced in conjunction with the butter churn.

The cream setting pan and dairy thermometer in combination will produce cream from 4 litres of milk. Depending on the fat content of the milk, approximately 0.5 litres of cream can be obtained from this and the minimum churn charge of 1 litre can be quickly separated. This can be used if no automatic cream separator such as FT15 (Disc Bowl Centrifuge) is available. Approximately 0.5kg of butter can be made from 1.5 litres of cream in the churn. During the churning process it is necessary to drain of the buttermilk and washing water and a sieve is provided to strain out any butter particles to be returned to the churn. Scotch Hands are required to 'work' the butter into a solid mass after churning. Water exuded during the working process is removed by the use of muslin cloth.

In order to carry out these processes, the following accessories are supplied:

- 1 x 4 litre capacity cream setting pan
- 1 x dairy thermometer
- 1 x nylon sieve
- 1 x pair of Scotch Hands
- 1 x roll butter muslin

Armfield Limited Bridge House West Street Ringwood Hampshire England BH24 1DY Tel: +44 (0)1425 478781 Fax: +44 (0)1425 470916 E mail: sales@armfield.co.uk URL: http://www.armfield.co.uk

USA Office:

Armfield Inc. 436 West Commodore Blvd (#2) Jackson NJ 08527 Tel: (732) 928-3332 Fax: (732) 928-3542 E mail: info@armfieldinc.com